

MODIFICATIONS

(Made at time the shoe is manufactured)

Duplicate cast (fiberglass negative returned)	
Make positive cast from old shoe	
Casting prosthetic foot	
Reconstruct broken casts	
Metatarsal pads (per foot cast modification)	
Pockets on insoles (per foot cast modification)	
Combination removable insole (dual density) consisting of 1/4" white 1/4 pink plastazote®	
Third layer, each additional layer, or substitution of material (other than 1/4" pink + 1/4" white plastazote®)	
Leather cover on insoles	
Filler for toe amputation (or for cosmetic purposes)	
Molded forefoot Amputation Filler	
Specialty leather uppers (nubuck, heavy duty or deerskin)	
Full or partial. Plastazote® or leather lining	
Straps with buckles (substituted for laces)	
Additional Velcro® straps	
Zipppers	
Speed laces (metal loops)	
Leather loop pull strap	
Bellowed Tongues	
Additional laces	
T-Strap (lateral or medial)	
Extended (double) medial or lateral counters	
Steel or fiberglass toes	
Caliper plate (your, reinforced and installed)	
Solid Stirrup (yours, reinforced and installed)	
Graphite Mortons extension plate or Graphite plate	
Steel shank (full sole bar)	
Welting (cemented)	
Heel elevation (crepe)	
Full elevation (crepe)	
Full elevation (cork, leather covered)	
Heel Flare (lateral or medial)	
Sole Flare (lateral or medial)	
Heel Wedge (lateral or medial)	
Sole Wedge (lateral or medial)	
Wide base	
Heavy duty ribbed sole	
Neoprene sole (oil resistant)	
Topy sole	
Leather sole	
Leather sole with outside heel	
Leather toe tips	
Toe and heel saver	
S.A.C.H. Heel	
Metatarsal bars	
Base excavation	
Water protector treatment (no guarantee)	
Rush service (shipped 3 working days from date order confirm)	

REPAIRS AND MODIFICATIONS

(Can be made to any existing molded shoe)

Full Repair: Balance, resole with standard 18 iron soling	
Heel Elevation Crepe includes new sole	
Elevation Platform Crepe includes new sole	
Elevation Platform Cork leather covered	
Rocker Soles	
Steel Shank includes new sole	
Reinforced Counter includes new sole	
Wedge or Flare new sole	
Rolled Padded Collar	
Replace Velcro	
Convert existing Velcro closure to lace	
Convert existing lace closure to Velcro strap	
Insoles from positive cast sent in	
Met Pads new top cover	
Toe filler added to existing insole	
New heels only	
New toe tips only	
Caliper Plate yours reinforced and installed includes new sole	

WE CAN DUPLICATE ANY STYLE SHOE

WE ACCEPT DIFFERENT KINDS OF CAST:

- PLASTER
- FIBERGLASS WRAP
- STS CASTING SOCK
- (PLEASE MAKE SURE THE TOES HAVE BEEN ACURATELY CAPTURED)

THE FOLLOWING INFORMATION WILL HELP US MAKE THE BEST SHOES POSSIBLE:

- WEIGHT BEARING TRACING
- MANUAL FOOT MEASUREMENTS (BALL & INSTEP)
- FLOOR REACTIONS IMPRINTS
- TRACING OF THE CURRENT INSERTS

12 DAYS TURNAROUND